

			Patient Information				
			Patient Full Legal Name (Last, First)				
1111 6th Avenue Des Moines, IA 50314 Phone 515-247-4471 Fax 515-643-8832			Date of Birth: / /		Sex: ☐ Male ☐ Female		
			Phone #		SSN #		
Location		Location ID	Address				
ADDRESS			City, State, Zip Code				
CITY		STATE	ZIP	Collection Date:			
PHONE			Collection Time:				
Ordering Provider (NPI):			Billing Information				
Ordering Provider (Print):			☐ Bill Account				
☐ Fax To # ☐ Call To #							
For MercyOne DSM Lab Use Only			Additional Tests:				
Date/Time Rcvd:		Initials:	RQ#:				
MRN #		Courier Tracking #					
Employee/Occupational Health Requisition Form							
	Allergen IgE Latex	LAB1737	Gold Clot/Spin/Refrigerate		Hepatitis B Surface Antibody	LAB472	Gold Clot/Spin/ Refrigerate
	Alpha Fetoprotein Tumor Marker	LAB559	Gold Clot/Spin/Aliquot/Freeze		Hepatitis B Surface Antigen	LAB8951	Gold Clot/Spin/Refrigerate
	Basic Metabolic Panel	LAB15	Green Spin/Refrigerate		Hepatitis C Antibody Screen	LAB10241	Gold Clot/Spin/Refrigerate
	Blood Typing (ABO/Rh)	LAB895	Lavender (Do Not Spin) Refrigerate		HIV Type 1 and 2	LAB10235	Gold Clot/Spin/Refrigerate
	C Reactive Protein HS	LAB150	Green Spin/Refrigerate		Lead Level, Venous	LAB98	Dark Blue (Do Not Spin) Refrigerate *Add'l Form Needed
	CA 125	LAB155	Gold ClotSpin/Refrigerate		Lipid Panel	LAB18	Green Spin/Refrigerate
	CA 15 3	LAB776	Red Clot/Spin/Aliquot/Refrigerate		Magnesium Level	LAB103	Green Spin/Refrigerate
	CA 19 9	LAB777	Gold Clot/Spin/Refrigerate		Mumps Virus IgG Antibody	LAB160	Gold Clot/Spin/Refrigerate
	Cadmium Level, Whole Blood	LAB4237	Dark Blue (Do Not Spin) Refrigerate		Neuron-Specific Enolase	LAB4623	Red Clot/Spin/Aliquot/Refrigerate
	Carcinoembryonic Antigen	LAB57	Gold Clot/Spin/Refrigerate		Potassium Level	LAB114	Green Refrigerate
	CBC Hemogram	LAB294	Lavender (Do Not Spin) Refrigerate		Prostate Specific Angiten (PSA) Diagnostic	LAB4010	Gold Clot/Spin/Refrigerate
	Cholesterol Total	LAB60	Green Spin/Refrigerate		Prostate Specific Antigen (PSA) Screen	LAB10243	Gold Clot/Spin/Refrigerate
	Chromium Creatinine Ratio Urine, Random	Misc.	Urine (Mayo: CRCRU) Refrigerate		PSA Total & Free	LAB171	Gold Clot/Spin/Aliquot/Freeze
	Chromium Level	LAB967	Dark Blue (Do Not Spin) Refrigerate		Protine (PT/INR)	LAB320	Light Blue Spin/Aliquot/Freeze
	Complete Blood Count	LAB294	Lavender (Do Not Spin) Refrigerate		Partial Thromboplastin Time (PTT)	LAB325	Light Blue Spin/Aliquot/Freeze
	Complete Blood Count w/Differential	LAB293	Lavender (Do Not Spin) Refrigerate		Quantiferon TB Gold Plus	LAB9520	QuantIFERON-TB Gold 4 Tube Kit (Gray, Green, Yellow, Purple)
	Comprehensive Metabolic Panel	LAB17	Green Spin/Refrigerate		Rubella IgG Antibody	LAB496	Green Spin/Refrigerate
	Drug Screen Collect	LAB6442	Urine Refrigerate		Rubeola IgG Antibody	LAB657	Gold Clot/Spin/Refrigerate
	Glucose, Fasting	LAB81	Green Spin/Refrigerate		T4 Free	LAB127	Green Spin/Refrigerate
	Glucose, Random	LAB82	Green Spin/Refrigerate		Thyroid Stimulating Hormone (TSH)	LAB129	Green Spin/Refrigerate
	HCG Blood Qualitative	LAB144	Green Spin/Refrigerate		Varicella Zoster Virus IgG Antibody	LAB162	Gold Clot/Spin/Refrigerate
	Heavy Metals Scrn, Whole Blood	LAB832	Dark Blue (Do Not Spin) Refrigerate		Vitamin D	LAB535	Gold Clot/Spin/Refrigerate
	Hemoglobin	LAB291	Lavender (Do Not Spin) Refrigerate				
	Hemoglobin A1c	LAB90	Lavender (Do Not Spin) Refrigerate				
	Hepatitis B Surface Ab Qual and Quant	LAB1337	Gold Clot/Spin/Refrigerate				

*Clot: Allow blood to sit undisturbed for 15-20min until solid clot has formed

*Spin: Centrifuge until red cells have separated from Serum/Plasma

*Aliquot: Pour off Serum/Plasma into separate tube with patient label & write Primary Tube color on the lable (ie. Red, Gold, etc.)